



# Cardiac Markers Control

A Fundamental Cardiac Marker For Your Immunoassay System

A routine tri-level cardiac marker quality control for the use as an assayed control material in monitoring the precision of CKMB, troponin, myoglobin and C-Reactive Protein (CRP).

- Easy-to-use liquid human serum-based control
- Frozen to maintain a **two-year shelf life**
- **Stable for 15 days** at 2-8°C after thawing
- Low cTnI level challenges low end of range where qualified for cardiac risk stratification

**More Diagnostics' Cardiac Markers Control** has assayed values for: **Siemens Dimension® and Beckman Access®** and is also appropriate for other automated immunoassay systems that correlate with these methods.

## ORDERING INFORMATION

CAT#	PRODUCT DESCRIPTION	QUANTITY
150	Tri-Level Cardiac Markers Controls	6x3mL
150-1	Level 1 Cardiac Markers Controls	6x3mL
150-2	Level 2 Cardiac Markers Controls	6x3mL
150-3	Level 3 Cardiac Markers Controls	6x3mL

This control is included in More Diagnostics Quality Assurance Program (QAP) program

- an interlaboratory statistical program
- continuously monitor the performance of products
- compare to your laboratory peer groups

This free service enables you to spot trends and maintain high-quality lab results.



For further information, please contact More Diagnostics:

[Sales@MoreDiagnostics.com](mailto:Sales@MoreDiagnostics.com) | 1-800-758-0978 | 2020 11th St, Los Osos, CA 93402