

**MORE DIAGNOSTICS • CARDIAC MARKERS CONTROL #150 • VALUE SHEET • EXPECTED VALUES**

EXPECTED VALUES	UNITS	LEVEL ONE 150-1		LEVEL TWO 150-2		LEVEL THREE 150-3	
		<b>LOT</b> 61941	EXPECTED RANGE	<b>LOT</b> 61942	EXPECTED RANGE	<b>LOT</b> 61943	EXPECTED RANGE
<b>TROPONIN</b>							
DIMENSION	ng/mL	0.05	0.02-0.08	1.91	1.05-2.77	27.35	19.15-35.55
BECKMAN P/N A98264	ng/mL	0.054	0.027-0.081	1.033	0.517-1.549	15.812	8.697-22.927
BECKMAN P/N A98143	ng/ml	0.054	0.027-0.081	1.139	0.569-1.709	18.374	10.106-26.642
<b>MYOGLOBIN</b>							
DIMENSION	ng/mL	44	31-57	179	125-233	367	257-477
BECKMAN ACCESS SYSTEMS	ng/mL	33.28	23.30-43.26	117.48	82.24-152.72	226.29	158.40-294.18
<b>CKMB (MASS)</b>							
DIMENSION	ng/mL	2.0	1.4-2.6	23.9	16.7-31.1	66.0	46.2-85.8
BECKMAN ACCESS SYSTEMS	ng/mL	3.52	2.46-4.58	31.02	21.71-40.33	72.35	50.65-94.05
<b>CRP</b>							
DIMENSION CCRP	mg/L	3.19	0.00-6.38	11.44	8.41-14.47	18.51	13.73-23.29
BECKMAN DXC CRP	mg/dL	0.26	0.18-0.34	0.99	0.87-1.11	1.63	1.30-1.96

**STABILITY:** FROZEN: THIS PRODUCT IS STABLE UNTIL EXPIRATION DATE ON THE VIAL WHEN KEPT AT LESS THAN -26°C.  
**STABILITY:** THAWED: THIS PRODUCT IS STABLE FOR 15 DAYS WHEN STORED TIGHTLY CAPPED AT 2 TO 8°C.

**DO NOT REFREEZE AFTER THAWING.** Rev. No. Original 2016/11 PN 150 010



**CCYY-MM-DD 2018-07-12 MORE DIAGNOSTICS, INC., LOS OSOS, CA. 93402 (800) 758-0978**